

Complex environment engineering diffractometer (CEED)

Conceptual design update

Update of performance characteristics by MC simulations

J. Šaroun, P. Beran, P. Lukáš, J. Pilch, P. Strunz,
Nuclear Physics Institute ASCR, v.v.i., Řež, Czech Republic


P. Šittner
Institute of Physics ASCR, v.v.i., Prague

New features tested by MC simulations


- Updated neutron guides geometry
- Radial collimators for small gauge volumes
- Double chopper for constant $\Delta\lambda/\lambda$ resolution
- Double frame option
- New virtual samples for “engineering” materials: duplex steel, TiAl alloy

CEED design


a very schematic view ...



CEED design - primary beam



Double chopper for pulse shaping



Test setup:

2-frames

$$\lambda_1 = 1.8 \text{ \AA}$$


$$\lambda_2 = 3.6 \text{ \AA}$$

total range:

$$\Delta\lambda = 0.9 \text{ .. } 4.5 \text{ \AA}$$

similar e.g. to D33/ILL
but limited source pulse
duration = limited band
width

Double chopper



RESTRAX view of the simulation model

Chopper characteristics

outer radius: 350 mm

frequency: 140 Hz

window angles: ~ 54 and 72 deg

distances:

high resolution (HR)


0.2 m ... $\tau = 50 \mu\text{s}$ at $\lambda=1 \text{ \AA}$

low resolution (LR)

1.0 m ... $\tau = 250 \mu\text{s}$ at $\lambda=1 \text{ \AA}$

Two choppers operating in “blind” mode:
chopper 1 closes when chopper 2 opens - provides flat $\Delta\lambda/\lambda$ - resolution curve

Sample area at 90° configuration



sample

D=8 mm, H=16 mm cylinder

RESTRAX view of the simulation model

Incident beam

free path to the sample: 400 mm

vertical: radial collimator

acceptance angle 3 deg

number of slits 17

length: 400 mm

horizontal: slit, variable width

Scattered beam

free path to the sample: 400 mm

vertical: no collimator

acceptance angle 30 deg


horizontal: radial collimator

acceptance angle 20 deg

number of slits 81

length: 800 mm

Sample area at 170° configuration



Incident beam
the same as at 90°

Scattered beam

vertical: no collimator


acceptance angle 30 deg

horizontal: no collimator

acceptance angle 10 deg

RESTRAX view of the simulation model


Brilliance transfer to the sample no radial collimators



Brilliance transfer as a function of wavelength

integration volume: $dS = 5 \times 10 \text{ mm}^2$, $d\Omega = 0.3 \times 0.3 \text{ deg}^2$

Time resolution




Nearly constant $\Delta\lambda/\lambda$ resolution as produced by the double chopper for different distances, $dc=0.2$ and $dc=1.0$ m.


Beam width 2 cm => worse resolution at short wavelengths

Gauge volume, $2\theta = 90^\circ$

projection along incident beam



projection on horizontal plane



Gauge volume is defined by the **two radial collimators**

- vertical (before the sample)
- horizontal (after the sample)

and input slit


Integrated volume:

$\sim 25 \text{ mm}^3$



free space around the sample: $R=400 \text{ mm}$

Gauge volume, $2\theta = 170^\circ$

projection along incident beam




projection on horizontal plane




Brilliance transfer

including choppers and radial collimators



integration volume: $dS = 3 \times 3 \text{ mm}^2$, $d\Omega = 0.002 \times 0.02 \text{ sr}$


Isotropic scattering



Conversion to d_{hkl} scale is the same as for a powder sample (event based analysis)

Resolution $\Delta d/d$

Simulation with calibration sample (equal F^2 and equidistant d values)



Resolution curves for high and low resolution modes

NOTES:


Resolution includes the contribution of beam divergence, sample size etc. (dominant at high resolution)

d -range is reduced for the *low resolution* mode due to the penumbra effect


Simulated experiment 1: TiAl

high resolution mode (as described before)
two wavelength frames
detector range 80 - 100 deg

simulated count rate




normalized with Vanadium




Simulated experiment 1: TiAl

comparison of detector angles *zoom at low d*

detector range 80 - 100 deg




detector range 165 - 170 deg



Simulated experiment 2: duplex steel

comparison of low and high resolution modes *zoom at low d*



wavelength band 0.9 - 2.7 Å
detector range 80 - 100 deg

Summary

- Detailed model of CEED instrument concept is available
- Small gauge volume feasible with radial collimators at ~ 40 cm from the sample (minor loss of intensity)
- Flat $\Delta\lambda/\lambda$ resolution curve, tunable with the double chopper
- **Penumbra effect - significant reduction of usable bandwidth** (mainly at low resolution and high λ) => *should the instrument be longer?*
- double frame mode: possible, but not easy to configure, avoid cross-talk effects etc.
- tests with virtual samples - basis for intercomparisons and testing data analysis in future